



# YOUSEQ

## ANNOUNCING OUR FIRST CE-IVD MARKED KIT MARCH 2021



We have recently received the fantastic news that our COVID-19 2 Gene MultiPlex qRT-PCR is officially CE IVD compliant. CE Marking indicates that it complies with the European In-Vitro Diagnostic Devices Directive (98/79/EC) and that the device may be legally commercialised in the EU.

Designed as a single tube 3 target format, our COVID-19 qPCR kit offers multiplex detection of two SARS-CoV-2 targets (E and RdRp genes plus an RNA control) for maximum diagnostic certainty. It's

compatible with all known mutants/variants and is ready to use with no resuspension steps required.

We are immensely proud of the team that have worked tirelessly over the last 6 months to achieve our COVID-19 qPCR CE Mark: Having our first CE-IVD Marked product is yet another step in our ongoing commitment to ensuring optimal product quality and delighting our customers. You can view our COVID-19 qPCR kit [Here](#)