

INTRODUCING OUR qPCR CUSTOM PANELS TO BUILD YOUR LAB DEVELOPED TEST

13th October 2022



YouSeq can now supply your laboratory with convenient, qPCR-reaction-ready, primers/probes dried in to 96 or 384 well plates.

We wanted to offer our customers a first-class product, designed using our proprietary AI driven methods, that significantly increase the efficiency of LDT development. Our new custom panels are incredibly convenient – simply add qPCR MasterMix and your nucleic acid sample for a perfect qPCR reaction.

All panels are designed to your own specification and validated by our expert team under ISO13485. Oligos can be supplied in a format most suited for each customer's instrumentation. There is also the flexibility to choose from our extensive library of targets or request custom development of a new plate.

All plates can be supplied alongside our proprietary MasterMix, Tetra™, which drives efficient amplification and probe cleavage of up to four targets at once. YouSeq can also provide controls separately or pre-dispensed into the oligo plate for your convenience. Furthermore, we now offer automated data analysis, powered by FastFinder for a complete high-throughput solution.

Please contact us to discuss your specific custom plate requirements, we would be delighted to help.

 FastFinder

