



YOUSEQ
DETECT MORE

MANUFACTURING SCIENTIST I

£22,500 - £23,500 Full Time

YouSeq is a biotechnology start-up. We make test kits that are used in hospitals and labs around the world to help fight COVID19, cancer and other diseases.

We are looking for a meticulous individual to join our team to help QC and manufacture qPCR and NGS test kits. The successful candidate will, after training, play a key role in our quality control process. On the surface our QC and manufacturing processes are simple: They all centre around putting chemicals in tubes. But to do it efficiently and at scale takes precision, organisation and brilliant record keeping.

YouSeq is a start-up and this role comes with all the fun that that brings. Your work-life will be varied and the company is growing fast. You'll need to be ready to pivot and move as the company grows.

The ideal candidate will be a meticulous perfectionist; good at organising processes and keeping accurate records.

Essential skills/experience:

- Use of pipettes to dispense precise volumes
- Demonstrable precision and organisation skills

Desirable skills:

- Knowledge of qPCR as a technology
- Experience of using qPCR as a technology
- Experience of working in ISO13485 manufacturing environment

Qualifications:

BSc or equivalent

Benefits:

- Casual dress
- Company pension
- Private medical insurance
- Buy and sell annual leave scheme
- On-site parking

If you would like to apply, please forward a covering letter and CV to hr@youseq.com