

STI qPCR PANELS

Powered by Tetra™ MasterMix

INTRODUCING OUR STI qPCR PANELS, DESIGNED WITH YOUR LAB IN MIND



YouSeq specialise in designing multiplex qPCR pathogen detection kits. We are now offering our customers the freedom and convenience to select exactly which of these targets they would like pre-dispensed into 96 well or 384 well plates to create larger panels.

You can choose from our current library of targets, or we can custom develop them for you. We can then design the plate layout for you to fit in the maximum number of samples to drive efficiency and detect more.

TARGET PATHOGENS AVAILABLE

Sexually Transmitted Infections

Atopobium vaginae	Herpes Simplex Virus 2
BVAB-2	Human Herpes Virus 6
Candida albicans	Mobiluncus curtisii
Candida parapsilosis	Mobiluncus mulieris
Candida tropicalis	Mycoplasma genitalium
Chlamydia trachomatis	Mycoplasma hominis
Cytomegalovirus	Neisseria gonorrhoeae
Enterococcus faecalis	Staphylococcus aureus
Escherichia coli	Streptococcus agalactiae (Grp B)
Gardnerella vaginalis	Streptococcus pyogenes (Grp A)
Haemophilus ducreyi	Treponema pallidum
High Risk HPV 16	Trichomonas vaginalis
High Risk HPV 18	Ureaplasma parvum
High Risk HPV 45	Ureaplasma urealyticum
Herpes Simplex Virus 1	Varicella zoster virus

