

QC LABORATORY TECHNICIAN

6-month contract

YouSeq is a biotechnology company that makes molecular testing kits that are used in hospitals and labs around the world to help fight cancer and other diseases.

This role is responsible for performing quality control and manufacturing processes to ensure the accurate and timely production of positive control templates and singleplex units. You will maintain a well organised laboratory environment by overseeing the stock take and maintenance of consumables, as well as managing the ordering and goods-in processes for QC lab items. You will also collaborate closely with the dispatch team to ensure the correct and timely supply of positive control templates for customers.

Job responsibilities include but are not limited to:

- Stock take and maintenance of positive control template and QC lab consumables
- Rack, lid and label components manufactured in the QC Labs
- Ordering and goods-in of QC lab related items
- Preparing gBlocks for dispatch
- Routine QC lab coat and mop head cleaning

Essential qualifications/skills:

- Excellent eye for detail

Desirable:

- Confidence working to quality management systems, e.g. ISO 9001/ISO 13485

Job type:

- Full time, initially a 6-month contract with the potential to become a permanent role

Work location:

- Laboratory based (Winchester)

Benefits:

- Casual dress
- Company events
- On-site parking

If you would like to apply, please forward a covering letter and CV to hr@youseq.com